



## DAD Software

*Beckman Instruments, Inc.*

ArrayView, a software option for the model 168 diode array detector is an easy-to-use, intuitive package that is accessed via a menu item in Gold software and reads Gold data files directly. It generates a 4-in-1 display of an isoabsorbance map, a spectrum, a 3-D plot, and a chromatogram. Each of these four windows can be expanded for closer examination. The interactive isoabsorbance map allows the user to click on any wavelength to immediately display the chromatogram at the wavelength selected and also the spectrum at the chosen time. These data can be saved, transferred into Gold software, and integrated as any other data file.