

Only Smart Diode Array Detection points out the impurities. In real time.

Smart Diode Array Detection shows you which peaks are impure. And how many components exist in each. Automatically, and in *real-time*.

Unique Purity Lines below each peak make it obvious. So you don't have to guess like you did with vague, unclear 3-D plots and isograms from older generation detectors. Or wait for time-consuming post-run analyses.

The next generation of diode array detection is here today. It's the System Gold™ 168 Diode Array Detector. It's Smart. And it's only from Beckman.

For a free electronic brochure, contact your Beckman Representative. Offices in major cities worldwide.

Call 800/742-2345 in the US.

Or write Beckman Instruments, Inc., Bioanalytical Systems Group, 4550 Norris Canyon Road, P.O. Box 5101, San Ramon, CA 94583.



BECKMAN

Australia, Sydney Austria, Vienna Canada, Toronto
France, Gagny Germany, Munich Hong Kong,
Aberdeen Italy, Milano Japan, Tokyo Mexico,
Mexico City Netherlands, Mijdrecht
Puerto Rico, Carolina Singapore,
South Africa, Johannesburg Spain,
Madrid Sweden, Bromma
Switzerland, Nyon Taiwan,
Taipei United Kingdom, High
Wycombe
©1990 Beckman Instruments,
Inc. AX90-1595B